


Supplementary data


Supplementary table 1 | **Fluorescent filters available in the Fluidigm system.**

Excitation filter	Emission filter
475-41 nm	525-25 nm
530-16 nm	570-30 nm
575-31 nm	630-30 nm
632-27 nm*	700-30 nm*


*optional filter





Supplementary figure 1 | **Reaction scheme of fluorescent assay for hydrogen peroxide detection on example of lactate oxidase.** Lactate is oxidised to pyruvate by lactate oxidase (LOX), generating hydrogen peroxide which reacts with a non-fluorescent probe AmplifluRed in the presence of horseradish peroxidase (HRP) to produce the fluorescent oxidation product, resorufin.



Supplementary figure 2 | **Resorufin fluorescence in the Fluidigm system and microplate reader.** Represented are ranges for linear responses captured in the two systems, subsequently used as resorufin standard curves.




Supplementary figure 3 | **Influence of buffer composition on fluorescent readout in the Fluidigm system.** (A) Fluorescent readouts for samples in buffers with and without addition of 0.01mg/ml BSA (final concentration), in presence of 5mM lactate. (B) Fluorescent readouts for samples in buffers with 0.01% and 0.1% concentration of detergents (Tween20 and Triton-X-100). Lactate oxidase (LOX) was used in these experiments. Error bars represent standard deviation of the data obtained with four replicates.




Supplementary figure 4 | **Distribution of fluorescent signal of the loading control fluorescein.** Signals from samples containing detergents Triton-X-100 and Tween20 in the concentration of 0.1% (left) and 0.01% (right).


A



B



C



Supplementary figure 5 | **Comparison of kinetic values obtained in the Fluidigm system and a plate reader.** Michaelis-Menten curves of three repeats in each system for the three tested oxidases: (A) Lactate oxidase. (B) Glutamate oxidase. (C) Glucose oxidase.